



Product Information Sheet

A1375 *Arabidopsis* Germination Medium

Synonym: *Arabidopsis* seed germination medium

Properties:

Form: Fine to Fluffy Powder
Appearance: White to Yellow Powder
Application: *Arabidopsis* seed germination
Solubility: Water
Typical Working Concentration: 12.67 g/L
Storage Temp: 2-6°C
Storage Temp of Stock Solution: Preparation of concentrated solutions is not recommended as insoluble precipitates may form.
Other Notes: Contains the macro- and micronutrients as described by Murashige and Skoog (1962).
pH = 3.5-4.5

Formula [mg/L]:

Ammonium Nitrate	825
Potassium Nitrate	950
Boric Acid	3.1
Cobalt Chloride 6H ₂ O	0.0125
Cupric Sulfate 5H ₂ O	0.0125
Na ₂ EDTA 2H ₂ O	18.63
Ferrous Sulfate 7H ₂ O	13.9
Magnesium Sulfate Anhydrous	90.35
Manganese Sulfate H ₂ O	8.45
Sodium Molybdate 2H ₂ O	0.125
Potassium Iodide	0.415

Potassium Phosphate Monobasic	85.0
Zinc Sulfate 7H ₂ O	4.3
Calcium Chloride Anhydrous	166.1
MES	500
Sucrose	10,000

Application Notes:

Arabidopsis thaliana (Col-0) and (Ler) possess a relatively low level of seed dormancy, and it is most often broken by stratification (2-6°C) in the dark for 2-4 days.

This formulation contains 1/2X MS (M524) with 0.5 g/L MES, and 10 g/L of sucrose and meets what the ABRC (Arabidopsis Biological Resource Center) recommends (<https://abrc.osu.edu/seed-handling>) for *Arabidopsis* seed germination. This formulation also contains enough nitrogen (20 mM in nitrate) and carbon (30 mM in sucrose) for proper seedling growth and development (Martin *et al.* 2002). As *Arabidopsis thaliana* is a model flowering plant, this type of formulation has been used often for germinating seeds (Denness *et al.* 2011; Stacey *et al.* 2006; Sorek *et al.* 2014).

If you prefer 1X MS for *Arabidopsis* germination, we would suggest using M530 (contains M519 + 0.5 g/L MES) and add 10 g/L of sucrose (Valvekens *et al.* 1988; Millet *et al.* 2010). To add Gamborg's vitamins add 1 mL/L of G219.

PhytoTechnology Laboratories®

P.O. Box 12205; Shawnee Mission, KS 66282-2205
Phone: 1-888-749-8682 or 1-913-341-5343; Fax: 1-888-449-8682 or 1-913-341-5442
Web Site: www.phytotechlab.com © 2014 PhytoTechnology Laboratories®



Product Information Sheet

References:

Denness L, McKenna JF, Segonzac C, Wormit A, Madhou P, Bennett M, Mansfield J, Zipfel C, and Hamann T. (2011) Cell Wall Damage-Induced Lignin Biosynthesis Is Regulated by a Reactive Oxygen Species- and Jasmonic Acid-Dependent Process in Arabidopsis. *Plant Physiol.* Jul 2011; 156(3): 1364–1374.

Martin T, Oswald O, and Graham IA (2002) Arabidopsis Seedling Growth, Storage Lipid Mobilization, and Photosynthetic Gene Expression Are Regulated by Carbon:Nitrogen Availability. *Plant Physiol.* Vol. 128 (2) pg. 472-481

Stacey MG, Osawa H, Patel A, Gassmann W, and Stacey G (2006) Expression analyses of Arabidopsis oligopeptide transporters during seed germination, vegetative growth and reproduction. *Planta* Vol. 223(2) pg. 291-305

Sorek N, Sorek H, Kijac A, Szemenyei HJ, Bauer S, Hématy K, Wemmer DE, and Somerville CR (2014) The Arabidopsis COBRA protein facilitates cellulose crystallization at the plasma membrane. *JBC* First published online 20Oct2014 doi: 10.1074/jbc.M114.607192

Valvekens D, Van Montagu M, and Van Lijsebettens M. (1988) Agrobacterium tumefaciens- mediated transformation of Arabidopsis thaliana root explants by using kanamycin selection. *PNAS* 85, 5536-5540

Millet YA, Danna CH, Clay NK, Songnuan W, Simon MD, Werck-Reichhart D, and Ausubel FM (2010) Innate Immune Responses Activated in Arabidopsis Roots by Microbe-Associated Molecular Patterns. *Plant Cell* Vol. 22 pg. 973-990