




Product Information Sheet

D259
DL-Dithiothreitol



Synonym: 1,4-Dithiothreitol
DTT; Cleland's Reagent
CAS: 3483-12-3
Formula: C₄H₁₀O₂S₂
Molecular Wt: 154.2

Properties

Form: Powder
Appearance: White Powder
Application: Molecular Biology
Solubility: Water
Typical Working
Concentration:
Storage Temp: 2 to 6° C
Storage Temp of
Stock Solution:
Other Notes:

Application Notes

References

Merck 13, 3412

Revised 3/2007

PhytoTechnology Laboratories, LLC™

P.O. Box 12205 • Shawnee Mission, KS • 66282-2205

Phone: 1-913-341-5343 or 1-888-749-8682 (U.S. Only) Fax: 1-913-341-5442 or 1-888-449-8682 (U.S. Only)

Web Site: www.phytotechlab.com

© 2007 PhytoTechnology Laboratories, LLC™