



## Product Information Sheet

### E330 Eriksson Vitamin Solution (1000x)

#### Properties

Form:	Liquid
Appearance:	Clear, Colorless Solution
Application:	Plant Tissue Culture
Solubility:	Miscible with Water
Typical Working Concentration:	1.0 mL/L of Medium
Storage Temp:	2 – 6° C
Storage Temp of Stock Solution:	N/A
Other Notes:	Contains the vitamins as described by Eriksson (1965). pH = 3.75 – 4.75

#### Formula (mg/L)

Glycine (Free Base)	2000
Nicotinic Acid (Free Acid)	500
Pyridoxine•HCl	500
Thiamine•HCl	500

#### Application Notes

Plant Tissue Culture Tested

#### References

Eriksson, T. 1965. *Physiol. Plant* 18:976-993.

Revised 2/2007

### PhytoTechnology Laboratories, LLC

P.O. Box 12205; Shawnee Mission, KS 66282-2205

Phone: 1-888-749-8682 or 913-341-5343; Fax: 1-888-449-8682 or 913-341-5442

Web Site: [www.phytotechlab.com](http://www.phytotechlab.com)

© 2007 PhytoTechnology Laboratories, LLC