

PhytoTechnology Laboratories, LLC™

"Helping To Build A Better Tomorrow Through Plant Science"™

Product Information Sheet

A126 Antioxidant Mixture

Properties

Form:	Fine to Fluffy Powder
Appearance:	White
Application:	Antioxidant
Solubility:	Water
Typical Working	0.25 g/L
Concentration:	0.25 g/L
Storage Temp:	Room Temperature
Storage Temp of	2-6° C
Stock Solution:	2-0 C
Other Notes:	Plant Tissue Culture Tested

Formula (mg/L)

L-Ascorbic Acid	100
Citric Acid	150

Application Notes

This product typically is made up as a stock solution and an explant is dipped into it before transfer to fresh medium.

References

Revised 3/2007

PhytoTechnology Laboratories, LLC™

P.O. Box 12205 • Shawnee Mission, KS • 66282-2205 Phone: 1-913-341-5343 or 1-888-749-8682 (U.S. Only) Fax: 1-913-341-5442 or 1-888-449-8682 (U.S. Only) Web Site: www.phytotechlab.com © 2007 *Phyto* Technology Laboratories, LLC[™]