



Product Information Sheet

M498 Malt Extract Agar

Properties

Form:	Powder
Appearance:	Off-White to Cream Powder
Application:	Bacterial Screening
Solubility:	Hot Water
Typical Working Concentration:	48.0 g/L
Storage Temp:	Room Temperature
Other Notes:	This product should be suspended in water and heated to boiling before dispensing and autoclaving.

Formula (mg/L)

Agar	15,000
Malt Extract	30,000
Peptone from Meat	3000

Application Notes

Typically used for the cultivation of yeasts and molds.

Revised Oct 2014

PhytoTechnology Laboratories, LLC

P.O. Box 12205; Shawnee Mission, KS 66282-2205

Phone: 1-888-749-8682 or 913-341-5343; Fax: 1-888-449-8682 or 913-341-5442

Web Site: www.phytotechlab.com

© 2007 PhytoTechnology Laboratories, LLC