



Product Information Sheet

O156 Orchid Maintenance Medium Contains Banana and Charcoal, Without Agar

Properties

Form:	Powder
Appearance:	Gray powder
Application:	Orchid Culture
Solubility:	Partially Soluble in Water
Typical Working Concentration:	57.31 g/L
Storage Temp:	2 – 6° C
Storage Temp of Stock Solution:	Preparation of concentrated solutions is not recommended as insoluble precipitates may form.
Other Notes:	pH = 4.5 – 5.5

Formula

(mg/L)

Ammonium Nitrate	825
Boric Acid	3.1
Calcium Chloride, Anhydrous	166
Cobalt Chloride•6H ₂ O	0.0125
Cupric Sulfate•5H ₂ O	0.0125
Na ₂ EDTA•2H ₂ O	37.3
Ferrous Sulfate•7H ₂ O	27.85
Magnesium Sulfate, Anhydrous	90.35
Manganese Sulfate•H ₂ O	8.45
Molybdic Acid (Sodium Salt) •2H ₂ O	0.125
Potassium Iodide	0.415
Potassium Nitrate	950

Potassium Phosphate, Monobasic	85
Zinc Sulfate•7H ₂ O	5.3
Activated Charcoal	2000
Banana Powder	30000
MES (Free Acid)	1000
myo-Inositol	100
Nicotinic Acid (Free Acid)	1
Peptone from Meat	2000
Pyridoxine•HCl	1
Sucrose	20000
Thiamine•HCl	10

Application Notes

Plant Tissue Culture Tested

Plant species: Epiphytic orchid species

This medium is a modification of our Orchid Maintenance/Replate Medium (P748). It lacks a gelling agent to permit its use in a liquid culture system. This medium contains banana powder to promote growth and rooting. The pH of this medium should be adjusted after adding the gelling agent and prior to sterilization

Revised 12/2012

PhytoTechnology Laboratories®

P.O. Box 12205; Shawnee Mission, KS 66282-2205

Phone: 888-749-8682 or 913-341-5343; Fax: 888-449-8682 or 1-913-341-5442

Web Site: www.phytotechlab.com

© 2012 PhytoTechnology Laboratories®