



Product Information Sheet

M5632 Deconstructed MS Group II (KNO₃) 50x Murashige & Skoog (MS)

Properties

Form:	Liquid
Appearance:	Colorless, Clear Liquid
Application:	Plant Tissue Culture
Solubility:	Miscible with Water
Typical Working Concentration:	20 mL/L
Storage Temp:	2-6°C
Other Notes:	Contains the potassium and nitrogen components as described by Murashige and Skoog (1962). pH = 4.5-7.0

Formula (mg/L)

Potassium Nitrate	95,000
-------------------	--------

Application Notes

Plant Tissue Culture Tested

This product contains only potassium and nitrogen macronutrients, other nutrients supplied via M5631, M5633, M5634, and M5635.

References

Murashige, T and F Skoog. 1962. A revised medium for rapid growth and bioassays with tobacco tissue cultures. *Physiol. Plant.* 15: 473-497.

Reed BM, Wada S, DeNoma J, Niedz RP (2013) Mineral nutrition influences physiological responses of pear *in vitro*. *In Vitro Cell.Dev.Biol.-Plant.* DOI 10.1007/s11627-013-9556-2