



Product Information Sheet

O753

Orchid Multiplication Medium Contains Agar, Without Charcoal Mother Flasking Medium III

Properties

Form:	Powder
Appearance:	White to Yellow powder
Application:	Orchid Culture
Solubility:	Hot Water
Typical Working Concentration:	32.30 g/L
Storage Temp:	2 – 6° C
Storage Temp of Stock Solution:	Preparation of concentrated solutions is not recommended as insoluble precipitates may form.
Other Notes:	pH = 4.75 – 5.75

Formula

(mg/L)

Ammonium Nitrate	825
Boric Acid	3.1
Calcium Chloride, Anhydrous	166
Cobalt Chloride•6H ₂ O	0.0125
Cupric Sulfate•5H ₂ O	0.0125
Na ₂ EDTA•2H ₂ O	37.3
Ferrous Sulfate•7H ₂ O	27.85
Magnesium Sulfate, Anhydrous	90.35
Manganese Sulfate•H ₂ O	8.45
Molybdc Acid, Sodium Salt•2H ₂ O	0.125
Potassium Iodide	0.415
Potassium Nitrate	950

Potassium Phosphate, Monobasic	85
Zinc Sulfate•7H ₂ O	5.3
Agar	7000
6-Benzylaminopurine (BA)	2.0
MES (Free Acid)	1000
myo-Inositol	100
α-Naphthaleneacetic Acid	0.5
Nicotinic Acid (Free Acid)	0.5
Peptone from Meat	2000
Pyridoxine•HCl	0.5
Sucrose	20,000
Thiamine•HCl	1

Application Notes

Plant Tissue Culture Tested

Plant species: Many epiphytic orchid species

This medium is a modification of our Orchid Multiplication Medium (P723) by the deletion of charcoal. This medium should work well for the propagation of Phalaenopsis and Dendrobium stem props.

Revised 12/2012

PhytoTechnology Laboratories®

P.O. Box 12205; Shawnee Mission, KS 66282-2205

Phone: 1-888-749-8682 or 1-913-341-5343; Fax: 1-888-449-8682 or 1-913-341-5442

Web Site: www.phytotechlab.com

© 2012 PhytoTechnology Laboratories®