



## Product Information Sheet

# KCl

**P730**  
**Potassium Chloride Solution 3.0 M**

Synonym: Electrode Filling/ Soaking Solution  
CAS: 7447-40-7  
Formula: KCl  
Molecular Wt: 74.55

### Properties

Form: Liquid  
Appearance: Colorless, Clear Solution  
Application: Electrode Filling/ Soaking  
Solubility: Miscible with Water  
Typical Working Concentration:  
Storage Temp: Room Temperature  
Storage Temp of Stock Solution: N/A  
Other Notes:

### Formula (g/L)

Potassium Chloride	223.65
--------------------	--------

### Application Notes

### References

Merck 13, 7704

Revised 3/2007

**PhytoTechnology Laboratories, LLC™**

P.O. Box 12205 • Shawnee Mission, KS • 66282-2205  
Phone: 1-913-341-5343 or 1-888-749-8682 (U.S. Only) Fax: 1-913-341-5442 or 1-888-449-8682 (U.S. Only)  
Web Site: [www.phytotechlab.com](http://www.phytotechlab.com) © 2007 PhytoTechnology Laboratories, LLC™