




Product Information Sheet

P770
Pluronic® F-68



Synonym: Poloxamer 188
CAS: 9003-11-6
Formula: $(C_3H_6O \cdot C_2H_4O)_x$
Molecular Wt: 8400 Avg.

Properties

Form: Powder
Appearance: White Powder
Application: Surfactant
Solubility: Water
Typical Working
Concentration:
Storage Temp: Room Temperature
Storage Temp of
Stock Solution:
Other Notes: Plant Tissue Culture Tested

Application Notes

References

Merck 13, 7645

Revised 3/2007

PhytoTechnology Laboratories, LLC™

P.O. Box 12205 • Shawnee Mission, KS • 66282-2205

Phone: 1-913-341-5343 or 1-888-749-8682 (U.S. Only) Fax: 1-913-341-5442 or 1-888-449-8682 (U.S. Only)

Web Site: www.phytotechlab.com

© 2007 PhytoTechnology Laboratories, LLC™