


Product Information Sheet

S807
L-Serine



Synonym: (S)-2-Amino-3-hydroxypropanoic Acid
CAS: 56-45-1
Formula: C₃H₇O₃N
Molecular Wt: 105.09

Properties

Form: Powder
Appearance: White Crystalline
Application: Amino Acid
Solubility: Water
Typical Working
Concentration:
Storage Temp: Room Temperature
Storage Temp of
Stock Solution:
Other Notes: Plant Tissue Culture Tested
Derived from Bacterial Fermentation of a Non-animal Source

Application Notes

References

Merck 13, 8534

Revised 3/2007

PhytoTechnology Laboratories, LLC™

P.O. Box 12205 • Shawnee Mission, KS • 66282-2205

Phone: 1-913-341-5343 or 1-888-749-8682 (U.S. Only) Fax: 1-913-341-5442 or 1-888-449-8682 (U.S. Only)

Web Site: www.phytotechlab.com

© 2007 PhytoTechnology Laboratories, LLC™