




## Product Information Sheet



**S663**  
**TRISODIUM CITRATE**  
**DIHYDRATE**

CAS: 6132-04-3  
Formula:  $C_6H_5Na_3O_7 \cdot 2H_2O$   
Molecular Wt: 294.10

### Properties

Form: Granular Powder  
Appearance: White  
Application: Generally used as a buffer  
Solubility: Soluble in Water  
Typical Working Concentration: N/A  
Storage Temp: Room Temperature  
Storage Temp of Stock Solution: Room Temperature

### Application Notes

Commonly used for the preparation of molecular biology buffers.  
Trisodium citrate dihydrate is a tribasic salt of citric acid.

**PhytoTechnology Laboratories®**

P.O. Box 12205; Shawnee Mission, KS 66282-2205

Phone: 1-888-749-8682 or 1-913-341-5343; Fax: 1-888-449-8682 or 1-913-341-5442

Web Site: [www.phytotechlab.com](http://www.phytotechlab.com)

© 2010 PhytoTechnology Laboratories®