



Product Information Sheet

S7536 Sabouraud Dextrose Broth

Properties

Form:	Powder
Appearance:	Cream to Tan Powder
Application:	Microbiology & Phytopathology
Solubility:	Soluble in Water
Typical Working Concentration:	50.0 g/L
Storage Temp:	Room Temperature
Storage Temp of Stock Solution:	Preparation of concentrated solution is not recommended as insoluble precipitates may form.

Formula (mg/L)

D-Glucose, Anhydrous	40,000
Peptone, from Meat	10,000

Application Notes

Used in the isolation and cultivation of yeasts, fungi, and bacteria.

References

Difco Manual, 1984. Tenth Edition.