



Product Information Sheet

T754 Terrific Broth Solution w/ Glycerol (1x)

Properties

Form: Liquid
Appearance: Light Tan to Tan
Application: Molecular Biology
Solubility: Miscible with Water
Typical Working Concentration: 1x
Storage Temp: 2-6°C
Other Notes: Sterile Filtered

Formula (mg/L)

Glycerol	4.0 mL
Potassium Phosphate, Dibasic, Anhydrous	9400
Potassium Phosphate, Monobasic, Anhydrous	2200
Tryptone	11,800
Yeast Extract	23,600

Application Notes

Terrific Broth is a nutrient media for use with the culture *Escherichia coli* for microbiology applications.