



## Product Information Sheet

---

**B236**  
**Buffer Solution, pH Reference**  
**Standard, pH 7.0 ± 0.02 at 25° C**

### Properties

Form: Liquid  
Appearance: Clear, Yellow  
Application: pH Buffer Standard  
Solubility:  
Typical Working  
Concentration:  
Storage Temp: Room Temperature  
Storage Temp of  
Stock Solution:  
Other Notes:

### Application Notes

This is used for pH meter calibration

### References

Revised 3/2007

**PhytoTechnology Laboratories, LLC™**

P.O. Box 12205 • Shawnee Mission, KS • 66282-2205

Phone: 1-913-341-5343 or 1-888-749-8682 (U.S. Only) Fax: 1-913-341-5442 or 1-888-449-8682 (U.S. Only)

Web Site: [www.phytotechlab.com](http://www.phytotechlab.com)

© 2007 PhytoTechnology Laboratories, LLC™