



Product Information Sheet

B237

**Buffer Solution, pH Reference
Standard, pH 10.0 ± 0.02 at 25° C**

Properties

Form: Liquid
Appearance: Clear, Blue
Application: pH Buffer Standard
Solubility:
Typical Working
Concentration:
Storage Temp: Room Temperature
Storage Temp of
Stock Solution:
Other Notes:

Application Notes

This is used for pH meter calibration

References

Revised 3/2007

PhytoTechnology Laboratories, LLC™

P.O. Box 12205 • Shawnee Mission, KS • 66282-2205

Phone: 1-913-341-5343 or 1-888-749-8682 (U.S. Only) Fax: 1-913-341-5442 or 1-888-449-8682 (U.S. Only)

Web Site: www.phytotechlab.com

© 2007 PhytoTechnology Laboratories, LLC™