

# CERTIFICATE OF ANALYSIS

Product Description	POTASSIUM PHOSPHATE , MONOBASIC, ANHYDROUS ACS REAGENT MONOPOTASSIUM PHOSPHATE
Product Number	P846
Lot Number	HKA0846054
Storage Temperature	Room Temperature
Molecular Weight	136.09 g/mol
Formula	$\text{KH}_2\text{PO}_4$
CAS Number	7778-77-0

## Physiochemical Specifications:

TEST	SPECIFICATION	RESULTS
Solubility	Soluble in Water @ 10 mg/mL	Passes
pH	4.0-5.0	4.8
Physical Appearance Color* Texture	PM-1 to PM-4 Crystalline Powder	PM-1 White Crystalline Powder
Solution Appearance Clarity Color	Clear Colorless	Clear Colorless
Average Time to Dissolve	For Information Only (Approx. 3 min)	Passes
Insolubles	None	Passes
Purity	Minimum 99.0%	99.41%
Loss on Drying	For Information Only	0.04%

\* Product color based upon comparisons between sample and standardized color wheel (Benjamin Moore® Color Preview™).

## Biological Testing:

Test Concentration: 170 mg/L

TEST SPECIFICATION	PLANT CELL LINE	RESULTS
Supports and/or facilitates plant growth and/or shoot proliferation in two or more plant tissue cultured lines with no morphological aberrations to plants	Begonia	Passes
	Tobacco Callus	Passes

The material described in this certificate is synthetic. No animal- or plant-derived components were used in the manufacture of this product.

PhytoTech Labs Inc. provides the above information intended to be used only as a guide to the appropriate handling of this material by a properly trained person. PhytoTech Labs Inc. shall not be held liable for any damage resulting from handling or from contact with the above product. This product is intended for **LABORATORY USE ONLY**. Our products may NOT BE USED as drugs, cosmetics, agricultural or pesticidal products, food additives or as household chemicals.

Date of Release: 14-May-2020

Recommended Shelf Life Date: May-2028



David S. Hart  
Technical Services Manager

**PhytoTech Labs Inc.**  
14610 W 106<sup>TH</sup> St. Lenexa, KS 66215  
Phone: 1-888-749-8682 or 1-913-341-5343  
Fax: 1-888-449-8682 or 1-913-341-5442  
Phytotechlab.com