

CERTIFICATE OF ANALYSIS

Product Description ZEATIN RIBOSIDE, *trans* ISOMER
Product Number Z899
Lot Number HWY0899036
Storage Temperature -20°C
Molecular Weight 351.4 g/mol
Formula C₁₅H₂₁N₅O₅
CAS Number 6025-53-2

Physiochemical Specifications:

TEST	SPECIFICATION	RESULTS
Solubility	Soluble in 1 N KOH @ 1mg/mL	Passes
Physical Appearance Color*	PM-1 to PM-4, 2143-70 to 2148-70, or 2025-70 to 2029-70	PM-4 White
Texture	Fine to Granular Powder	Granular Powder
Solution Appearance Clarity Color	Clear Colorless to Slight Yellow Tint	Clear Colorless
Average Time to Dissolve	For Information Only (Approx. 1 min)	Passes
Insolubles	None	Passes
FTIR	Conforms to Structure	Passes
Purity by HPLC	Minimum 98.0%	98.28%
Loss on Drying	For Information Only	0.1%

* Product color based upon comparisons between sample and standardized color wheel (Benjamin Moore® Color Preview™).

Biological Testing:

Test Concentration: 2.0 mg/L

TEST SPECIFICATION	PLANT CELL LINE	RESULTS
Supports and/or facilitates plant growth and/or shoot proliferation in two or more plant tissue cultured lines with no morphological aberrations to plants	Achimenes	Passes
	Tobacco Callus	Passes

The material described in this certificate is synthetic. No animal- or plant-derived components were used in the manufacture of this product.

PhytoTech Labs Inc. provides the above information intended to be used only as a guide to the appropriate handling of this material by a properly trained person. PhytoTech Labs Inc. shall not be held liable for any damage resulting from handling or from contact with the above product. This product is intended for **LABORATORY USE ONLY**. Our products may NOT BE USED as drugs, cosmetics, agricultural or pesticidal products, food additives or as household chemicals.

Date of Release: 4-June-2021

Recommended Shelf Life Date: June 2024



David S. Hart
Technical Services Manager

PhytoTech Labs Inc.
14610 W 106TH St. Lenexa, KS 66215
Phone: 1-888-749-8682 or 1-913-341-5343
Fax: 1-888-449-8682 or 1-913-341-5442
Phytotechlab.com