




Product Information Sheet

D253

**Dichloroisocyanuric Acid,
Sodium Salt**



Synonym: NaDCC

CAS: 2893-78-9

Formula: C₃Cl₂N₃O₃Na

Molecular Wt: 219.9

Properties

Form: Powder

Appearance: White Granular

Application: Plant Tissue Disinfectant

Solubility: Water

Typical Working
Concentration: 2 to 5 grams per liter

Storage Temp: Room Temperature

Storage Temp of
Stock Solution: Room Temperature

Other Notes: Plant Tissue Culture Tested

Application Notes

This product typically is made up as a stock solution and plant tissue is soaked or rinsed in it before sterile water rinses and transfer to fresh medium. The amount of disinfecting time is dependant upon tissue type being used, and it is the responsibility of the purchaser to determine the appropriateness of this product for the specific plants that are being cultured and applications that are being used.

PhytoTechnology Laboratories, LLC™

P.O. Box 12205 • Shawnee Mission, KS • 66282-2205

Phone: 1-913-341-5343 or 1-888-749-8682 (U.S. Only) Fax: 1-913-341-5442 or 1-888-449-8682 (U.S. Only)

Web Site: www.phytotechlab.com

© 2014 PhytoTechnology Laboratories, LLC™