



## Product Information Sheet

# L5240 LB Broth, Miller Modification, ADP-Free

### Properties

Form:	Powder
Appearance:	Cream to Tan Powder
Application:	Bacterial & Fungal Screening; Microbiology & Phytopathology
Solubility:	Soluble in Water
Typical Working Concentration:	29 g/L
Storage Temp:	2 - 6° C
Storage Temp of Stock Solution:	Preparation of concentrated stock solution not recommended as insoluble precipitates may form.
Other Notes:	This medium is free of animal derived products (ADP).

Formula	(mg/L)
Peptone, from Soy	14,000
Sodium Chloride	10,000
Yeast Extract	5,000

### Application Notes

This medium is used to grow *Escherichia coli* for molecular biology applications. This is a rich medium that promotes rapid growth of pure cultures of recombinant strains.

### References

- Lennox, ES. 1955. Transduction of linked genetic characters of the host by bacteriophage P1. *Virology*. 1:190-206.
- Luria, SE, and JW Burrous. 1955. Hybridization between *Escherichia coli* and *Shigella*. *J. Bacteriol.* 74:461-476.
- Luria, SE, JN Adams, and RC Ting. 1960. Transduction of lactose-utilizing ability among strain of *E. coli* and *S. dysenteriae* and the properties of the transducing phage particles. *Virology*. 12:348-390.
- Miller, JH 1972. *Experiment in molecular genetics*. Cold Spring Harbor Laboratory, Cold Spring Harbor, New York.
- Sambrook, J, EF Fritsch, and T Maniatis. 1989. *Molecular Cloning: A Laboratory Manual*, 2nd Edition. Cold Spring Harbor Laboratory, Cold Spring Harbor, New York.

Rev. 7/2009

### PhytoTechnology Laboratories®

P.O. Box 12205; Shawnee Mission, KS 66282-2205

Phone: 1-888-749-8682 or 1-913-341-5343; Fax: 1-888-449-8682 or 1-913-341-5442

Web Site: [www.phytotechlab.com](http://www.phytotechlab.com)

© 2009 PhytoTechnology Laboratories®