




Product Information Sheet



M5561
Manganese Chloride Tetrahydrate, USP

CAS: 13446-34-9
Formula: $\text{MnCl}_2 \cdot 4\text{H}_2\text{O}$
Molecular Wt: 197.91

Properties

Form: Powder
Appearance: Pink to Red Crystalline
Application: Plant Micronutrient
Solubility: Water
Typical Working Concentration: See Application Note
Storage Temp: Room Temperature
Storage Temp of Stock Solution: Refrigerate to retard growth of microorganisms in non-sterile solutions.
Other Notes: Plant Tissue Culture Tested

Application Notes

This is the form of manganese that is contained in Hoagland's Modified Basal Salts (Prod. No. H353) at 1.81 mg/L.

References

Merck 13, 5751

Revised 5/2010

PhytoTechnology Laboratories®

P.O. Box 12205 • Shawnee Mission, KS • 66282-2205

Phone: 1-913-341-5343 or 1-888-749-8682 (U.S. Only) Fax: 1-913-341-5442 or 1-888-449-8682 (U.S. Only)

Web Site: www.phytotechlab.com

© 2010 PhytoTechnology Laboratories®