



Product Information Sheet

N5890 **Nutrient Agar, ADP-Free**

Properties

Form:	Powder
Appearance:	Cream to Tan Powder
Application:	Bacterial & Fungal Screening; Microbiology & Phytopathology
Solubility:	Partially Soluble in Water
Typical Working Concentration:	23.8 g/L
Storage Temp:	Room Temperature
Storage Temp of Stock Solution:	N/A
Other Notes:	Contains no animal-derived products (ADP) Preparation of concentrated solution is not recommended as insoluble precipitates may form.

Formula (mg/L)

Agar	15,000
Peptone, from Soy	8,800

Application Notes

Typically used in the cultivation of various bacteria, fungi, and yeasts.

Rev. 7/2009

PhytoTechnology Laboratories®

P.O. Box 12205; Shawnee Mission, KS 66282-2205

Phone: 1-888-749-8682 or 1-913-341-5343; Fax: 1-888-449-8682 or 1-913-341-5442

Web Site: www.phytotechlab.com

© 2009 PhytoTechnology Laboratories®