



Product Information Sheet

HNO₃

N617
Nitric Acid Solution 1.0 N

Synonym: Aquafortis
CAS: 7697-37-2
Formula: HNO₃
Molecular Wt: 63.01

Properties

Form: Liquid
Appearance: Colorless, Clear Liquid
Application: Solvent; pH Adjustment
Solubility: Miscible with Water
Typical Working
Concentration:
Storage Temp: Room Temperature
Storage Temp of
Stock Solution: N/A
Other Notes: Plant Tissue Culture Tested

Application Notes

References

Merck 13, 6608

Revised 3/2007

PhytoTechnology Laboratories, LLC™

P.O. Box 12205 • Shawnee Mission, KS • 66282-2205

Phone: 1-913-341-5343 or 1-888-749-8682 (U.S. Only) Fax: 1-913-341-5442 or 1-888-449-8682 (U.S. Only)

Web Site: www.phytotechlab.com

© 2007 PhytoTechnology Laboratories, LLC™