



Product Information Sheet

P6780 Potassium Sodium Tartrate, Tetrahydrate

Synonyms: Rochelle Salt, Seignette Salt, L(+)-Tartaric acid
potassium Sodium salt
CAS: 6381-59-5
Formula: $C_4H_4KNaO_6 \cdot 4H_2O$
Mol. Weight: 282.22

Properties

Form: Powder
Appearance: White to Off-White
Solubility: Soluble in Water
Application: Seed Testing, Molecular Biology
Storage Temp: Room Temperature
Typical Working
Concentration: 34.6g/100 mL
Other Notes: Fehling's Reagent (B)

Application Notes

A component of Fehling's solution, used to differentiate between water-soluble carbohydrate and ketone functional groups, and as a test for monosaccharides. Fehling's solution can be used to determine whether a carbonyl-containing compound is an aldehyde or a ketone.

Fehling's solution, is composed of 2 separately prepared solutions, added together in equal parts.

Fehling Soln A consists of 7 g of hydrated copper (II) sulfate (Prod# C375) dissolved in 100 mL of DI water.

Fehling Soln B is made by dissolving 34.6 g of potassium sodium tartrate and 24 g of sodium hydroxide (Prod # P672) in 100 mL of DI water. This solution is referred to as a general test for monosaccharides (Nigam, 2007) and should always be made fresh.

References

Nigam A. (2007) Lab Manual of Biochemistry. Tata McGraw-Hill Education, ISBN: 0070617678. Pg. 21-22.

Unrau AM (1964) Carbohydrates of Centrosema Seed. The Constitution of an Arabino Galactan. *Canadian Journal of Chemistry*, 42. Pp. 916-925.

PhytoTechnology Laboratories®

P.O. Box 12205; Shawnee Mission, KS 66282-2205

Phone: 1-888-749-8682 or 1-913-341-5343; Fax: 1-888-449-8682 or 1-913-341-5442

Web Site: www.phytotechlab.com

© 2012 PhytoTechnology Laboratories®