



Product Information Sheet

P732 Peptone Water, Buffered (1x)

Properties

Form: Liquid
Appearance: Clear, Light Tan to Tan
Application: Molecular Biology
Solubility: Miscible with Water
Typical Working Concentration: 1x
Storage Temp: 2-6°C
Other Notes: Buffered to An approximate pH of 7.2, Sterile Filtered

Formula (mg/L)

Peptone from Meat	10,000
Sodium Chloride	5000
Sodium Phosphate, Dibasic	3500
Potassium Phosphate, Monobasic	1500

Application Notes