




Product Information Sheet

P738

Polyoxyethylenesorbitan Monooleate



Synonym: Tween® 80; Polysorbate 80

CAS: 9005-65-6

Formula: $\text{C}_{64}\text{H}_{124}\text{O}_{26}$

Molecular Wt: 1309.63

Properties

Form: Liquid

Appearance: Clear, Yellow to Amber, Viscous Liquid

Application: Surfactant

Solubility: Miscible with Water

Typical Working

Concentration:

Storage Temp: Room Temperature

Storage Temp of

Stock Solution: N/A

Other Notes: Plant Tissue Culture Tested

Application Notes

References

Merck **13**, 7664

Revised 3/2007

PhytoTechnology Laboratories, LLC™

P.O. Box 12205 • Shawnee Mission, KS • 66282-2205

Phone: 1-913-341-5343 or 1-888-749-8682 (U.S. Only) Fax: 1-913-341-5442 or 1-888-449-8682 (U.S. Only)

Web Site: www.phytotechlab.com

© 2007 PhytoTechnology Laboratories, LLC™