



Product Information Sheet

P772 Potato Dextrose Agar

Properties

Form:	Powder
Appearance:	Off-White to Cream Powder
Application:	Culture Microorganisms
Solubility:	Hot Water
Typical Working Concentration:	39.0 g/L
Storage Temp:	Room Temperature
Other Notes:	This product should be suspended in water and heated to a boil if it is to be dispensed prior to autoclaving.

Formula (mg/L)

Agar	15,000
D-Glucose, Anhydrous	20,000
Potato Powder	4000

Application Notes

References

PhytoTechnology Laboratories, LLC

P.O. Box 12205; Shawnee Mission, KS 66282-2205
Phone: 1-888-749-8682 or 913-341-5343; Fax: 1-888-449-8682 or 913-341-5442
Web Site: www.phytotechlab.com © 2007 PhytoTechnology Laboratories, LLC